



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

ACETAMINOPHEN

Test Code: ACTMN

Specimen Type: Serum

Tube Type / Collection Container: BB (Brick color) or Red Top

Specimen Volume: 1 ML

Specimen Handling: Centrifuge and pipette serum into plastic screw cap transport tube.

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Serum must be removed from cells and pipette into plastic transport tube.

Department: Chemistry

