



One Hospital Drive • Lewisburg, PA 17837  
Laboratory Services - (570)522-2510

## ACTH

**Test Code:** ACTH

**Specimen Type:** Plasma

**Tube Type / Collection Container:** Lavender K3 EDTA, Prechilled

**Specimen Volume:** 0.1 ML

**Specimen Handling:** Centrifuge and pipette plasma into plastic screw cap transport tube.

**Specimen Storage:** Frozen

**Specimen Collection Notes:** Collect full ice chilled Lav tube; immediately after venipuncture place in ice bath, centrifuge asap, transfer plasma into a plastic transport tube, freeze immediately

**Department:** Reference Laboratory

**Reference Lab Code:** 004440



Collect full ice chilled Lav tube; immediately after venipuncture place in ice bath, centrifuge asap, transfer plasma into a plastic transport tube, freeze immediately

