

One Hospital Drive • Lewisburg, PA 17837 Laboratory Services - (570)522-2510

AFP Tumor Marker

Test Code:

Specimen Type: Serum

Tube Type / Collection Container: SST (Serum Separator Tube)

Specimen Volume: 1 ML

Specimen Handling: Centrifuge

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: This Test is for MALES and NON-

PREGNANT Females ONLY. Not intended for serial monitoring; One

time determination only.

Department: Reference Laboratory

Reference Lab Code: 002253

