



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

ALBUMIN, FLUID

Test Code: ALBFL

Specimen Type: Fluid

Tube Type / Collection Container: Dark Green; Sodium Heparin

Specimen Volume: 1 ML

Specimen Handling: Place fluid into dark green tube immediately after collection.

Specimen Storage: Refrigerate 2-8 C Stability 48hrs

Specimen Collection Notes: Fluid/Specify site-source

Department: Chemistry

