

One Hospital Drive • Lewisburg, PA 17837 Laboratory Services - (570)522-2510

ALDOSTERONE, serum

Test Code: ALDOS

Specimen Type: Serum

Tube Type / Collection Container: SST (Serum Separator Tube)

Specimen Volume: 1 ML

Specimen Handling: Centrifuge

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Only with physician approval the patient should have no isotopes administered for 24hrs prior to venipuncture, should be taken off medications for at least 3 weeks prior to sample collection and should be on a normal sodium diet. Patient should be ambulatory for a least 30min. before blood collection. / State patient position-supine/upright at time of collection.

Department: Reference Laboratory

Reference Lab Code: 004374

