



One Hospital Drive • Lewisburg, PA 17837  
Laboratory Services - (570)522-2510

## **Aldosterone, urine**

**Test Code:** MISC

**Specimen Type:** Urine

**Tube Type / Collection Container:** 24 hr Urine Container

**Specimen Volume:** Submit Entire 24 HR Collection

**Specimen Handling:** 1ML Aliquot

**Specimen Storage:** Refrigerate 2-8 C

**Specimen Collection Notes:** Only with physician approval; Patient should be on a diet containing 3g sodium per day for at least 2 weeks and preferably 30 days prior to testing and no isotopes administered 24hrs prior to venipuncture. / Container: Plastic 24hr w/1 g boric acid perserv. Refer to 24 hour urine collection procedure in this manual

**Department:** Reference Laboratory

**Reference Lab Code:** 004291



**Preservative:** 1 g Boric Acid