



One Hospital Drive • Lewisburg, PA 17837  
Laboratory Services - (570)522-2510

## AMIODARONE

**Test Code:** AMIO

**Specimen Type:** Serum

**Tube Type / Collection Container:** Brick Color Top OR Red Top

**Specimen Volume:** 2 ML

**Specimen Handling:** Centrifuge and pipette serum into plastic screw cap transport tube.

**Specimen Storage:** Refrigerate 2-8 C

**Specimen Collection Notes:** SST Tubes will be rejected. Test includes:  
Noramiodarone Metabolite

**Department:** Reference Laboratory

**Reference Lab Code:** 706705



OR

