



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

ANTI-Xa FACTOR

Test Code: AXA

Specimen Type: Plasma

Tube Type / Collection Container: Light Blue; 3.2% Buffered Sodium Citrate

Specimen Volume: 2.7 ML Draw to fill line on tube-Short draws will be rejected.

Specimen Handling: Centrifuge at 3000rpm within 1hr of collection and pipette plasma into plastic screw cap transport tube.

Specimen Storage: Ambient / Room Temp. Stability 2hrs at Room Temp.

Department: Chemistry

