



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

B. Pertussis (Smear and Culture)

Test Code: MCB

Specimen Type: Nasal

Tube Type / Collection Container: Nasal Swab / Slides

Specimen Volume: Swab AND 3-4 Slides


Specimen Handling: AirDry, Don't Fix

Specimen Storage: Ambient / Room Temp.

Specimen Collection Notes: 3-4 Microscope slides of material from nasal swab. Air dry. Do NOT Fix.

Department: Microbiology / Reference Laboratory

Reference Lab Code: 188425



3-4 Air Dried Slides DO NOT FIX