



One Hospital Drive • Lewisburg, PA 17837  
Laboratory Services - (570)522-2510

## BRONCHIAL BRUSH CULTURE

**Test Code:** PROCAT

**Specimen Type:** Bronchial Brush

**Tube Type / Collection Container:** 1 ML Tryp. Soy Broth

**Specimen Handling:** Place Bronchial Brush into Trypticase Soy Broth

**Specimen Storage:** Ambient / Room Temp.

**Specimen Collection Notes:** Transport ASAP.

**Department:** Microbiology

Bronchial Brush is placed into TSB



Brush should be sent to the Lab  
Immediately if site does not have  
Trypticase Soy Broth.



Place brush in TSB if  
available at collection site.