



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

CBC W/ AUTO DIFF

Test Code: CBC

Specimen Type: Whole Blood

Tube Type / Collection Container: Lavender K3 EDTA

Specimen Volume: Optimum 5 ML; Minimum 2 ML

Specimen Handling:

Specimen Storage: Ambient / Room Temp.

Specimen Collection Notes: Prefer 1 FULL lav top tube, but a minimum of 1/3rd full may be acceptable. Do not centrifuge/Whole Blood stable 24hrs Room Temp (can run Retic count off of same tube)

Department: Hematology

