



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

CPK, TOTAL

Test Code: CK

Specimen Type: Serum

Tube Type / Collection Container: SST (Serum Separator Tube)

Specimen Volume: 1 ML

Specimen Handling: Centrifuge

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Must draw discard tube prior to collecting a CPK specimen/Write draw off on appropriate tube and submit with CPK tube/Serum stable 48hrs refrigerated

Department: Chemistry

