



One Hospital Drive • Lewisburg, PA 17837  
Laboratory Services - (570)522-2510

## CULTURE-URINE

**Test Code:** URINE

**Specimen Type:** Urine

**Tube Type / Collection Container:** Sterile Container

**Specimen Volume:** Special Urine CULTURE transfer tube (Vaccuette)

**Specimen Handling:** Vaccuette Urine tube contains preservative.

**Specimen Storage:** Refrigerate 2-8 C

**Specimen Collection Notes:** Prefer 1st morning void, refrigerated no more than 8 hrs. Stable 24 hrs. Transfer a portion of the urine specimen from collection container to a CULTURE transport tube using a urine transfer kit. Submit both specimens. Specify Clean Catch or Catheterization collection. Contact lab for kit.

**Department:** Microbiology



**Collect Sample in Sterile cup (Clean Catch or Cath) and Transfer into “Vacuette Urine Tube” for transport to Laboratory.**