



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

CULTURE-UROGENITAL

Test Code: UROGE

Specimen Type: Swab

Tube Type / Collection Container: Swab with Stewarts media or eSwab

Specimen Handling: Collect with appropriate swab transport ASAP

Specimen Storage: Ambient / Room Temp.

Specimen Collection Notes: Transport As Soon As Possible

Department: Microbiology



Used for Routine (non anaerobic)
Urogenital culture.

OR



Used when BOTH Routine
AND Anaerobic Urogenital
cultures are requested.