



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

CULTURE-WOUND

Test Code: WOUND

Specimen Type:

Tube Type / Collection Container: Swab with Stewarts media or eSwab

Specimen Handling: Collect with appropriate swab transport ASAP

Specimen Storage: Ambient / Room Temp.

Specimen Collection Notes: Transport ASAP.

Department: Microbiology



Used for Routine (non anaerobic)
Urogenital culture.

OR

Used when BOTH Routine
AND Anaerobic Urogenital
cultures are requested.

