



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

CULTURE-YERSINIA

Test Code: YER

Specimen Type: FECES (STOOL)

Tube Type / Collection Container: Stool Collection Kit

Specimen Volume: 1-2g/ ML

Specimen Handling: Use Stool Collection Kit. Follow the enclosed patient instructions.

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: 1-2g (thumbnail size portion) Stool or Rectal swab submitted in a C & S (Clean vial-White pictured below) or Cary Blair transport vial (Orange pictured below). Refrigerate immediately after collection.

Department: Microbiology

