



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

CALCIUM, IONIZED

Test Code: CAION

Specimen Type: Serum

Tube Type / Collection Container: SST (Serum Separator Tube);
Unopened

Specimen Volume: 1 ML

Specimen Handling: Centrifuge; **Do Not Open tube.**

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Place a piece of tape over the stopper with
“Ionized Calcium – Do Not Open” written on it.

Department: Reference Laboratory

Reference Lab Code: 82330



Use Paper Tape Over Top of Tube Stopper