



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

COMPLEMENT C-3

Test Code: COMP3

Specimen Type: Serum

Tube Type / Collection Container: SST (Serum Separator Tube)

Specimen Volume: 1 ML

Specimen Handling: Centrifuge

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Refrigerate; Allow to clot at Room Temp. 15 to 30 mins.

Department: Reference Laboratory

Reference Lab Code: 006452

