



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

COMPLEMENT TOTAL (CH50)

Test Code: CH50

Specimen Type: Serum

Tube Type / Collection Container: SST (Serum Separator Tube)

Specimen Volume: 1 ML

Specimen Handling: Centrifuge

Specimen Storage: Frozen

Specimen Collection Notes: Allow to clot at room temp for 15-30 mins. Remove serum after centrifugation, place in transport tube and **MUST** freeze immediately

Department: Reference Laboratory

Reference Lab Code: 001941

