

One Hospital Drive • Lewisburg, PA 17837 Laboratory Services - (570)522-2510

CORTISOL URINE, FREE by HPLC

Test Code: CORTU

Specimen Type: Urine

Tube Type / Collection Container: 24 hr Urine Container

Specimen Volume: 100 ML Aliquot

Specimen Handling: Submit Entire 24 HR Collection

Specimen Storage: Room Temperature

Specimen Collection Notes: 1 gram Boric Acid or HCL; 10 ML minimum Refer to 24 hour urine collection procedure in this manual; Refrigerate during Collection

Department: Reference Laboratory

Reference Lab Code: 004432

2	
1	A manage straight
ł	and the second s
	The second second
	MARKING AND INT
	Berginskinger at
	Themail and
	ATTREAMEN (P

Preservative: 1 gram Boric Acid or HCL **No Preservative is also Acceptable**