

One Hospital Drive • Lewisburg, PA 17837 Laboratory Services - (570)522-2510

CREATININE CLEARANCE

Test Code: CRECL

Specimen Type: Urine & Serum

Tube Type / Collection Container: 24 hr Urine Container and SST (Serum

Separator Tube)

Specimen Volume: 10 ML

Specimen Handling: Centrifuge SST; Submit entire 24hr urine collection

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: 24hr Urine & Serum & patient's height & weight required/Urine must be refrigerated/SST stable 48hrs refrigerated Refer to 24 hour urine collection procedure in this manual

Department: Chemistry





No Preservative Needed