



One Hospital Drive • Lewisburg, PA 17837  
Laboratory Services - (570)522-2510

## CULTURE-CATHETER TIP

**Test Code:** CATHTIP

**Specimen Type:** Catheter Tip

**Tube Type / Collection Container:** Sterile Container

**Specimen Storage:** Ambient / Room Temp.

**Specimen Collection Notes:** Aseptically place catheter tip into a sterile container. Transport ASAP. Stable 24 hours at room temperature

**Department:** Microbiology



Actual Cup Appearance May Vary – This Specimen Requires Use of a Sterile Cup