



One Hospital Drive • Lewisburg, PA 17837  
Laboratory Services - (570)522-2510

## CULTURE-GC

**Test Code:** GC

**Specimen Type:** Vaginal or Urethral Swab

**Tube Type / Collection Container:** Thayer Martin Jembec or Swab

**Specimen Handling:** BBL Culture Swab Plus onto Jembec

**Specimen Storage:** Ambient / Room Temp.


**Specimen Collection Notes:** Transport ASAP.

**Department:** Microbiology

Male collection Swab is shown  
on Top, Female Below.



After swab collection, gently roll  
/ rub swab onto "Jembec Media"  
activate pouch as indicated  
(Instructions provided).



Jembec Media