



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

DISOPYRAMIDE (NORPACE)

Test Code: DISOP

Specimen Type: Serum

Tube Type / Collection Container: Brick Color Top Or Red Top

Specimen Volume: 1 ML

Specimen Handling: Centrifuge and pipette serum into plastic screw cap transport tube.

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: NO SST / Centrifuge & remove from cells within 2 hrs of collection / use plastic transport tube

Department: Reference Laboratory

Reference Lab Code: 007864



OR

