



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

DRUG SCREEN, 10 Panel

Test Code: TDS

Specimen Type: Urine

Tube Type / Collection Container: Clean Container

Specimen Volume: 10 ML

Specimen Handling:

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Urine stable 48 hours refrigerated

Department: Urinalysis



Actual specimen cup may appear different – This test does not require a sterile cup.