



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

FIBRIN SPLIT PRODUCTS

Test Code: FSP

Specimen Type: Serum

Tube Type / Collection Container: Special – Dk. Blue Hemagaurd closure tube with a Yellow label. (See Below)

Specimen Volume: 3 ML

Specimen Handling: Tube must set 30 – 60 mins. room temp before centrifugation.

Specimen Storage: Maintained at Room Temp while testing-No Add on testing can be requested.

Specimen Collection Notes: Must be collect at ECH hospital in special tube. All other tubes types are rejected.

Department: Chemistry

