

One Hospital Drive • Lewisburg, PA 17837 Laboratory Services - (570)522-2510

HEAVY METALS PROFILE I - BLOOD Lead, Mercury and Arsenic

Test Code: HEAV

Specimen Type: Whole Blood

Tube Type / Collection Container: Dark Blue Trace Metal Tube EDTA

Specimen Volume: 7 ML

Specimen Handling:

Specimen Storage: Ambient / Room Temp.

Specimen Collection Notes: Must submit original tube; Maintain at room temp.; Clotted spec.rejected; Patient should avoid seafood & red wine 72 hrs prior to collection.

Department: Reference Laboratory

Reference Lab Code: 042580

