



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

HOMOCYSTEINE

Test Code: HOMO

Specimen Type: Serum

Tube Type / Collection Container: Lavender, Dk Green or SST (Serum Separator Tube)

Specimen Volume: 2 ML

Specimen Handling: Centrifuge

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: 12hr fasting preferred/Specimen must be spun, after fully clotted if SST is used (within 30 mins of collection) and then immediately refrigerated

Department: Reference Laboratory

Reference Lab Code: 706994

