



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

KOH PREPARATION

Test Code: KOH

Specimen Type: Vaginal Swab, Nail, or Skin Scraping

Tube Type / Collection Container: Brick top/Red top or Clean Container or Swab in saline.

Specimen Volume:

Specimen Handling:

Specimen Storage: Ambient / Room Temp.

Specimen Collection Notes:

Department: Microbiology



in



or



or



Swab desired area and place into saline added to one of above tubes OR Place clippings (nails/skin into clean cup. Sterile Cup NOT Required.