



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

KLEIHAUER-BETKE

Test Code: BB

Specimen Type: Whole Blood

Tube Type / Collection Container: Lavender K2 EDTA

Specimen Volume: 3 ML

Specimen Handling:

Specimen Storage: Ambient / Room Temp.

Specimen Collection Notes: Do not centrifuge/Whole blood specimen must be maintained at room temperature and arrive in ECH same day as collected.

Department: Blood Bank

