



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

Latex IGE

Test Code: MISC

Specimen Type: Serum

Tube Type / Collection Container: SST (Serum Separator Tube)

Specimen Volume: 1 ML

Specimen Handling: Centrifuge

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Drawn tube must contain a minimum volume of 8.5ML prior to centrifugation. Specimen will be rejected if hemolyzed or grossly lipemic.

Department: Reference Laboratory

Reference Lab Code: 650390



