



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

LEGIONELLA PNEUM, ANTIBODY PANEL G,M

Test Code: LEGIO

Specimen Type: Serum

Tube Type / Collection Container: 2 SST (Serum Separator Tube)

Specimen Volume: 0.4 ML

Specimen Handling: Centrifuge

Specimen Storage: Ambient / Room Temp.

Specimen Collection Notes: Send paired serum samples, **one during the acute stage, and the second, three to six weeks later**, during the convalescent stage.

Department: Reference Laboratory

Reference Lab Code: 164616



Drawn at Acute Stage



Then

3 Weeks Later