



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

LIDOCAINE (Xylocaine®)

Test Code: LIDO

Specimen Type: Serum

Tube Type / Collection Container: Brick Color top, Red top; no additive
OR Dk Green

Specimen Volume: 1 ML

Specimen Handling: Centrifuge and pipette serum / Plasma into plastic
screw cap transport tube within 2HRS

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Draw 12Hrs after initiating therapy then every
24Hrs after.

Department: Reference Laboratory

Reference Lab Code: 007013



OR



OR



Transfer to

