

One Hospital Drive • Lewisburg, PA 17837 Laboratory Services - (570)522-2510

Manual Differential

Test Code: DIFFM

Specimen Type: Whole Blood

Tube Type / Collection Container: Lavender K3 EDTA

Specimen Volume: 2 ML

Specimen Handling:

Specimen Storage: Ambient / Room Temp.

Specimen Collection Notes: Prefer 1 FULL lav top tube, but a minimum of 1/3rd full may be acceptable. Do not centrifuge. Whole Blood stable 24hrs Room Temp

Department: Hematology

