



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

OVA AND PARASITE, STOOL

Test Code: OP

Specimen Type: Feces

Tube Type / Collection Container: Stool Collection Kit

Specimen Volume: 1-2g/ ML

Specimen Handling: Formalin & PVA Vials

Specimen Storage: Room Temperature

Specimen Collection Notes: 1-2g (thumbnail size portion) Stool submitted in the formalin & PVA vials. Refrigerate immediately after collection. 3-4mL Sputum or 15 mL Urine in transport is also acceptable.

Department: Microbiology / Reference Laboratory

Reference Lab Code: 008623

