



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

PERIPHERAL BLOOD SMEAR; REVIEWED BY PATHOLOGIST

Test Code: PERIP

Specimen Type: Whole Blood

Tube Type / Collection Container: Lavender K3 EDTA

Specimen Volume: 2 ML

Specimen Handling:

Specimen Storage: Ambient / Room Temp.

Specimen Collection Notes: If CBC performed at another facility: Prefer L top or 2 unstained smears accompanied by CBC report. Must include Diagnosis under consideration.

Department: Hematology



Pathologist Review