



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

PINWORM PREP

Test Code: PINWP

Specimen Type:

Tube Type / Collection Container: SWUBE Paddle

Specimen Volume:

Specimen Handling:

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Hold the paddle by the cap and remove it from the tube. Separate the buttocks and press the tacky surface against several areas of the perianal region. Replace the paddle in the tube for transport to the laboratory; specimens should be refrigerated if examination is to be delayed for more than one day.

Department: Microbiology

