

One Hospital Drive • Lewisburg, PA 17837 Laboratory Services - (570)522-2510

Rotovirus Detection by EIA

Test Code: MCB

Specimen Type: Feces

Tube Type / Collection Container: Clean Container

Specimen Volume: 2 ML liquid or semiformed stool

Specimen Handling:

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: 2ml liquid stool or 1g semiformed stool in screw cap container-3-5days after onset / Refrig immediately-To lab within 72 hrs.

Department: Microbiology / Reference Laboratory

Reference Lab Code: 006866



Actual Cup may vary in Appearance – Sterile Cup NOT Required.