



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

SALICYLATES

Test Code: SALIC

Specimen Type: Serum

Tube Type / Collection Container: Brick Color Top or Red Top

Specimen Volume: 2 ML

Specimen Handling: Centrifuge and pipette serum into plastic screw cap transport tube.

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Serum must be removed from cells within 2 hours of collection. Serum stable 48 hours refrigerated

Department: Chemistry



OR

