

One Hospital Drive • Lewisburg, PA 17837 Laboratory Services - (570)522-2510

SED RATE, AUTOMATED

Test Code: ESR

Specimen Type: whole Blood

Tube Type / Collection Container: ESR Tube

Specimen Volume: 1.2 ML

Specimen Handling: Mix well

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Stable 72 hours refrigerated. Fill to line; invert slowly 6-8 times for complete mixing.

Department: Hematology

