



One Hospital Drive • Lewisburg, PA 17837  
Laboratory Services - (570)522-2510

## **RAPID "A " STREP SCREEN, THROAT PCN SENSITIVE PATIENT**

**Test Code:** RSSAP

**Specimen Type:** Throat Swab

**Tube Type / Collection Container:** Swab

**Specimen Volume:**

**Specimen Handling:**

**Specimen Storage:** Refrigerate 2-8 C

**Specimen Collection Notes:** Order this test for patients that are sensitive to Pennicillin. Polyester or Rayon swab in Stuart media

**Department:** Hematology



Stuarts Media