



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

TROPONIN I

Test Code: TROPI

Specimen Type: Plasma

Tube Type / Collection Container: MINT GREEN- PST (Lithium Heparin) Plasma Separator Tube

Specimen Volume: 2 ML

Specimen Handling: Centrifuge

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Full Light Green w/Sodium Heparin. Must centrifuge within 6hrs of collection. Specimen stable 72 hours refig. If test is ordered for cardiac related symptoms, specimen should be transported to the lab same day or patient should be directed to the hospital lab for collection of specimen.

Department: Chemistry

