



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

ROUTINE URINALYSIS

Test Code: UA

Specimen Type: Urine

Tube Type / Collection Container: Clean Container / Sterile Container

Specimen Volume: 10 ML

Specimen Handling:

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Prefer 1st morning void, refrigerated no more than 8 hrs. Stable 24 hrs. Specify Clean Catch or Catheterization collection.

Department: Urinalysis



Appearance of Actual Cup May Vary- Use of Sterile Cup Recommended for Clean Catch and Catheterization Collections as Culture may be Indicated.