



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

URINALYSIS With CULTURE IF INDICATED

Test Code: UACI

Specimen Type: Urine

Tube Type / Collection Container: Sterile Container

Specimen Volume: 10 ML & filled transfer tube

Specimen Handling:

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Prefer 1st morning void. Stable 24 hours refrigerated. Transfer 1 ML of the urine specimen from collection container to a CULTURE transport tube using a urine transfer kit. Submit both specimens. Specify Clean Catch or Catheterization collection. Contact lab for kit.

Department: Urinalysis



Appearance of Actual Cup May Vary- USE OF STERILE CUP IS REQUIRED.