



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

VANCOMYCIN, PEAK

Test Code: VANCP

Specimen Type: Serum

Tube Type / Collection Container: Brick Color Top OR Red Top

Specimen Volume: 1 ML

Specimen Handling: Centrifuge and pipette serum into plastic screw cap transport tube.

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Serum must be removed from cells within 4 hours of collection. Specimen stable 48 hours refrigerated. Collect one hour after dose. Peak and Trough Vancomycin levels should be collected around the third dose.

Department: Chemistry



OR



