



One Hospital Drive • Lewisburg, PA 17837  
Laboratory Services - (570)522-2510

## VITAMIN C

**Test Code:** VITC

**Specimen Type:** Serum

**Tube Type / Collection Container:** SST (Serum Separator)

**Specimen Volume:** 2 ML

**Specimen Handling:** Centrifuge and protect from light.

**Specimen Storage:** Frozen

**Specimen Collection Notes:** Draw in chilled tube, protect from light, keep on ice, after clot, separate serum ASAP, pipette into transport tube, Prot. from Lght & Freeze immediately. Identify specimen w/patient name directly on the container AND on the outside of the aluminum foil.

**Department:** Reference Laboratory

**Reference Lab Code:** 001479



Draw in chilled tube-transfer serum to amber transport or wrap