

One Hospital Drive • Lewisburg, PA 17837 Laboratory Services - (570)522-2510

## VON WILLEBRAND PROFILE (Factor VIII); Includes VonWill ACT & Antigen

**Test Code:** MISC

**Specimen Type:** Plasma

Tube Type / Collection Container: 2 Light Blue; Buffered Sodium Citrate

3.2%(x2)

**Specimen Volume:** (2) 2 ML tubes (4 ML)

**Specimen Handling:** Centrifuge and pipette plasma into plastic screw cap transport tube.

or the cr

**Specimen Storage:** Frozen

**Specimen Collection Notes:** Call Lab for special collection instructions-Patient to avoid warfarin (Coumadin) therapy for 2 wks and heparin therapy for 2 days prior to test-Do Not stop taking med's without physician's approval-Do NOT draw from an arm with a heparin lock or heparinized catheter

**Department:** Reference Laboratory

Reference Lab Code: 084715

