



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

CLOSTRIDIUM DIFFICILE TOXIN A&B

Test Code: CLOSD

Specimen Type: Feces

Tube Type / Collection Container: Plain Vial

Specimen Volume: 1-2g/ ML

Specimen Handling: Plain Vial

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: 1-2g (thumbnail size portion) Unformed Stool submitted in a plain transport vial. Refrigerate immediately after collection. Performance characteristics not established for patients under 2yrs old.

Department: Microbiology



Clean Container-Sterile Cup is NOT Required.