



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

GONNORHEA(GC) BY DNA PROBE

Test Code: GCDNA

Specimen Type: Swab or Urine

Tube Type / Collection Container: Swab or Urine in Probetec Media

Specimen Volume: Swab in media or 30 ML urine

Specimen Handling: Request special kit for collections

Specimen Storage: Refrigerate 2-8 C

Specimen Collection Notes: Endocervical or urethral collection requires a special swab. Urine collection requires a sterile container, but DO NOT use the clean catch container or wipes. Urine should be a "dirty" first morning void. Please see GC and Chlamydia Urine Collection Protocol and GC and Chlamydia Swab Collection Protocol in the Customer Manual.

Department: Microbiology



Endocervical/Urethral :
Place "media swab"
directly into media, for
Urine sources: Void
into cup and transfer to
media as directed by
the GC/Chlamydia
Urine Collection
Protocol.

