



One Hospital Drive • Lewisburg, PA 17837
Laboratory Services - (570)522-2510

PHENYLALANINE QUANTITATIVE

Test Code: PHNYL

Specimen Type: Plasma, Frozen

Tube Type / Collection Container: Green top (Heparin) tube

Specimen Volume: 0.5 ML

Specimen Handling: Centrifuge and pipette serum into plastic screw cap transport tube, Freeze immediately.

Specimen Storage: Frozen

Specimen Collection Notes: Serum must be removed from cells immediately after collection and frozen. Specimens received thawed will be rejected.

Department: Reference Laboratory

